


AI/ML Services Hands On Lab:

In these labs, we will explore building and exploring scalable, serverless, machine learning applications using the AI/ML services. A key advantage of these services is the ability to build powerful ML application *without needing to be an expert machine learning practitioner*.

The following are services that we will explore from the console. You can click on the hyperlinks to go to the relevant AI/ML Service.

(**note** that in order for the links to the services to work, you must be logged into an AWS account)

1. Go to the Amazon [Comprehend Medical](#) console. Confirm that there is already input text, and then click `Analyze`. Explore the output and tabs that are present.
 - **Optionally**, enter your own text to see Comprehend Medical's capabilities on your text.
2. Go to the [Amazon Comprehend Console](#). Confirm that there is already input text, and then click `Analyze`. Explore the output and tabs that are present.
 - Within the input text, enter the text "My name is Josh. J as in Jungle, O as in Orange, S as in Sushi, H as in Horace.". Click on the `PII` section to see the results.
 - Within the input text, enter the text "I am exceptionally happy with my physician". Click on the Sentiment section to see the result.
3. Go to the [Amazon Textract Console](#). Confirm that there is already input document, and inspect the results.
 - Your Lab workshop comes with a file called `sample_patient_note.png`. Upload it Textract by clicking and dragging it to the console (it will take Textract a few seconds to analyze it). Inspect the results.
 - **Optionally** download the results by clicking `Download results`. Unzip the file, and inspect the JSON and other files created.
4. Go to the [Amazon Rekognition Console](#). Click on `Text in image`.
 - Your Lab workshop comes with a file called `chest-xray.jpg`. Upload it to the console to process it. Inspect the results.
5. We will now deploy a simple application that combines some of the AI/ML services we have learned about. The architecture of this application is:



Please follow these instructions to deploy and test the workshop:

- Go to the [AWS CloudFormation Console](#). Your Lab workshop comes with a file called ``cfn_workshop_ai_ml_services.yaml``. Click ``Upload a template file`` and upload this file.
- Enter any Stack name you wish; **do not change the Parameters**. Keep clicking next until you are asked to confirm the creation of the template. Click the checkbox, and then click ``Next``.
- The stack takes a few minutes to deploy. When it is done, go to ``Outputs``. It will show 2 buckets that are created. Note the name of your InputBucket (e.g. ``my-stack-notificationbucket-12iww3f3cdmf``) and your ResultsBucket (e.g. ``my-stack-outputbucket-wdxyeyczf7wa``).
- Go to the [Amazon S3 Console](#) find your InputBucket. Your Lab Workshop comes with a file called ``sample_patient_note.png``. Upload the document, wait a few seconds, and go to the ResultsBucket. Download and inspect the results.